



Release Notes

Floodlight MS SDK

Software version: 1.1.0 | Document version 1.0



Roche Molecular Systems, Inc.
2881 Scott Boulevard
Santa Clara, California 95050
United States of America

Date of Publication 28-04-2021

© 2021 F. Hoffmann-La Roche Ltd

All rights are reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without the prior written consent of the copyright owner.

Copyrights and all other proprietary rights in any software and related documentation ("Software") made available to you rest exclusively with Roche or its licensors. No title or ownership in the Software is conferred to you. Use of the software is subject to end user license agreement.

To the maximum extent permitted by law, you shall not decompile and/or reverse engineer the software or any part thereof.

Not for diagnostic, monitoring or therapeutic purposes or in any other manner for regular medical practice.

Preface

You use the Floodlight MS SDK to integrate the Floodlight MS assessments into third-party applications.

Unsupported Phones

The Floodlight MS SDK is not supported on the following phones:

Brand	Model Name
Alcatel	7
Alcatel	Tetra
Alcatel	ZIP LTE
Huawei	P20 Lite
Huawei	P Smart (2018)
Huawei	Mate 10 Lite
Huawei	Y6 (2018)
Huawei	P10
LG Electronics	Aristo 2
LG Electronics	K20 Plus
LG Electronics	Risio 3
LG Electronics	Tribute Dynasty
LG Electronics	Rebel 4
LG Electronics	Stylo 3 Plus
Motorola	Moto E5 Play
Motorola	Moto E4
Samsung	Galaxy J6

Known Limitations and Workarounds

The following issues are identified as known issues with workarounds in the current release of Floodlight MS SDK.

Floodlight MS Mobile Application crashes if you toggle the font

On Android phones, if the Floodlight MS Mobile Application is running in the background and you toggle the font, the app crashes. As a workaround, if you are using an Android device, do not use your device Settings to toggle the font while the Floodlight MS app is running.

Results are not getting displayed for Symptom Tracker

On iOS phones, if you use separate devices and perform the Symptom Tracker activity multiple times in a day, you will see incorrect results on your Dashboard and the Historical Result screen. You can perform the Symptom Tracker activity only once per day.

Daily Results are not displayed

If you cancel the Symptom Tracker activity before completing it, you will not see your results in the Daily Results screen. To view your results, complete the Symptom Tracker activity.

Data for the 365 day in the past is not visible in the PDF

The Historical Data exported as a PDF includes data for only 364 days. To view the data for the 365 days, use the Historical Data screen on your phone.

Test results are not displayed instantly in the bar charts

After you complete an activity, wait for a few seconds before you try to see the results of your activity in the Historical Results screen.

Inconsistent result status is displayed

If you manually change the date and time on your phone to a past date and time, you may not see your Historical Data and the past results on the Dashboard. To view your Historical Data and past results, reset the date and time to the Automatic setting on your phone.

Historical Data is not displayed correctly if the Date and Time on the phone is set to a date in the past

If you manually change the date and time on your phone to a past date and time, you may not see your Historical Data and the past results on the Dashboard. To view your Historical Data and past results, reset the date and time to the Automatic setting on your phone.

U Turn Test is not marked as done if you did not complete enough turns

If you completed the U-Turn activity but "Your Data" screen shows the message "This activity was not performed today," it indicates that you were not able to complete a minimum number of U-Turns in the given time. Please go through the user manual to understand the instructions and re-perform the activity.

PDF print-out displays data in UTC + 0 time zone

If you are exporting data to PDF to share it with your provider, the data will be in UTC. The data is also available in the charts screen and your provider also has access to it through the Floodlight MS Healthcare Professional's Portal.

Floodlight MS Mobile Application crashes if you change the device date and time from past to present

On an iOS phone, if you manually change the date and time on your phone to multiple past dates and times, you may not be able to access the Floodlight MS Mobile Application. To resolve this issue, reset the date and time to the Automatic setting on your phone.