



Release Notes

Floodlight MS Healthcare Professional's Portal

Software version: 1.1.0 | Document version 1.0



Roche Molecular Systems, Inc.
2881 Scott Boulevard
Santa Clara, California 95050
United States of America

Date of Publication 28-04-2021

© 2021 F. Hoffmann-La Roche Ltd

All rights are reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without the prior written consent of the copyright owner.

Copyrights and all other proprietary rights in any software and related documentation ("Software") made available to you rest exclusively with Roche or its licensors. No title or ownership in the Software is conferred to you. Use of the software is subject to end user license agreement.

To the maximum extent permitted by law, you shall not decompile and/or reverse engineer the software or any part thereof.

1 Preface

The Floodlight MS Healthcare Professional's Portal allows you to:

- View data collected by your patients using the Floodlight MS Mobile Application between clinic visits
- Review how MS impacts their cognition, hand function, and gait
- Improve the conversation around their care

1.1 Intended Use

The Floodlight Healthcare Professional's (HCP) Portal is intended to display the results of various assessments from the Floodlight Mobile Application per end-user and for several users over time. The HCP is solely responsible for interpreting the results.

1 Known Limitations and Workarounds

The following issues are identified as known issues with workarounds in the current release of Floodlight MS Healthcare Professional's Portal.

Y-axis overlays two bars with results when viewing patient chart with data for more than 90 days

If you have more than 90 days of patient data, the Y-axis may overlap some of the values on the far left. To resolve this issue, you can scroll the chart to see the data hidden by the Y-axis.

Automatic log-out doesn't work on some pages

When you first set up your account in the Floodlight MS Healthcare Professional's Portal, the automatic log out feature does not work. In this case, click Log out on the Navigation Panel at end of your first session.

Months on the chart are not translated when you switch the language on the patient test result (chart) view

If the provider changes the language when the test results are displayed in the Floodlight MS Healthcare Professional's Portal, the months in the chart are not translated. In this case, the provider can select the required language before navigating to the particular Patient's result page or refresh the page if changing the selected language while on the Patient's result page.