



Release Notes

Floodlight™ MS Healthcare Professional's Portal

Software version: 1.3 | Document version 5.0



Roche Molecular Systems, Inc.

2881 Scott Boulevard

Santa Clara, California 95050

United States of America



9342346001

Date of Publication 12-07-2021

Please contact Roche to request a printed version of this document.

© 2021 F. Hoffmann-La Roche Ltd

All rights are reserved. Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without the prior written consent of the copyright owner.




Copyrights and all other proprietary rights in any software and related documentation ("Software") made available to you rest exclusively with Roche or its licensors. No title or ownership in the Software is conferred to you. Use of the software is subject to end user license agreement.

To the maximum extent permitted by law, you shall not decompile and/or reverse engineer the software or any part thereof.

Not for diagnostic, monitoring or therapeutic purposes or in any other manner for regular medical practice.

Symbols used

The following table displays the symbols used in this Release Notes.

	Product manufacturer: The name and address of the application's manufacturer.
	Reference Number: Indicates the Reference or Catalog number for this user manual.
	Warnings and Precautions: Highlights information that is critical for optimal performance of the system. May also indicate that loss of data or invalid data could occur if the precautions or instructions are not observed.

Preface

The Floodlight™ MS Healthcare Professional's Portal allows you to:

- View data collected by your patients using the Floodlight™ MS Mobile Application between clinic visits
- Review how MS impacts their cognition, hand function, and gait
- Improve the conversation around their care

Intended Use

The Floodlight™ MS Healthcare Professional's Portal is intended to display the results of various assessments from the Floodlight™ MS Mobile Application per end-user and for several users over time. The HCP is solely responsible for interpreting the results.

Warnings and Precautions



Read the explanation for the data displayed in a patient's chart to understand the meaning of an increase and decrease of a metric value as well as the Y-axis labels. For more information, see the **Review a patient record** section in the Floodlight™ MS Healthcare Professional's Portal user manual.



Refer to the **Supported Platforms** (on page 7) section to ensure that the hardware and software requirements are met.

What's New?

This release of the Floodlight™ MS Healthcare Professional's Portal has been updated with new features to help enrich the conversation between you and your care team and offer a more convenient experience. Check it out!

- Improved look and feel and data export options for the Healthcare Providers.
- Notification for upcoming maintenance, recalls, and updates.
- Collect First Name, Last Name, and DOB for each patient at the time of sign up and include a patient identifier in the data export.
- Add the **Show Password** button to the **Log In** and **Change Password** screens to make it easier to debug log in issues.

Supported Platforms

The following table displays the supported platforms for the Floodlight™ MS Healthcare Professional's Portal.

Platforms	Version
DHP: dhp-api-gateway	4.9
DHP: dhp-gql-query-svc	4.2
floodlight-data-processor	4.3
DHP: dhp-patient-mgmt-svc	4.19
Floodlight™ MS SDK: iOS	1.3
Floodlight™ MS SDK Android	1.3
Floodlight™ MS SDK Backend	1.3

Known Limitations and Workarounds

The following issues are identified as known issues with workarounds in the current release of Floodlight™ MS Healthcare Professional's Portal.

Automatic log-out does not work on some pages

When you first set up your account in the Floodlight™ MS Healthcare Professional's Portal, the automatic log out feature does not work. In this case, click **Log out** on the Navigation Panel at the end of your first session.

Months on the charts are not translated when you switch the language on the patient test result (chart) view

If the provider changes the language when the test results are displayed in the Floodlight™ MS Healthcare Professional's Portal, the months in the chart are not translated. In this case, the provider can select the required language before navigating to the particular Patient's result page or refresh the page if changing the selected language while on the Patient's result page.

Team Owner can invite only registered users to their Team

In the Floodlight™ MS Healthcare Professional's Portal, a Team Owner can invite only a registered user as a Team Member. If you need to add a Team Member who is not already registered, please contact your system administrator or Roche. to first register the user.

Password field does not accept underscore (_) as a valid special character

In the Floodlight™ MS Healthcare Professional's Portal, the **Password** field does not accept underscore (_) as a valid special character. To avoid this error, instead on an underscore (_), use another special character in your password.

Known Limitations without present workarounds

The following limitations are unresolved in the current release of Floodlight™ MS Healthcare Professional's Portal. Resolutions or workarounds will be communicated in future release notes.

- The Team Member invitation email is sent in the language set for the sender instead of the language set for the recipient of the email.